

FOREVERMED®

Newly upgraded 90kv



Integrated design

Compact appearance
Light structure



Clear imaging

High efficiency



Small in size

Light in weight
Easy to operate



Low radiation level



Portable Digital X Ray High-frequency

YJ-BW905/BW9010

Application fields: Emergency room; Operating room; Sports medicine; Orthopedics; ICU; Nursing home; Physical examination; Ambulance; Rural medical visits; Field training; Animal husbandry; Pet hospital




TEL/WhatsApp: +86-15649625550
Email:ellen@yjcompany.com

Henan Forever Medical.Co.,Ltd.




Applicable to both humans and animals

To patient and operator with the lowest possible dose.



High frequency



Clarity



Efficient



Digital Flat Panel Detector
17×17-inch

Optional human-use or animal-use system

Compact design

Taking images from multiple angles

Optional features:

1. Exposure handbrake
2. Adjustable support bracket
3. Built-in computer system



Product model	YJ-BW905	YJ-BW9010
power	450W	900W
kV	40~90kV	40~90kV
mA	0.5~5mA	1~10mA
X-ray exposure time	0.1~2S	0.1~2S
Cooling time	Thirty times the exposure time	Sixty times the exposure time
focus size	0.8	0.8
Anode heat capacity	7000J	7000J
Total filtering	1.8mm Al	1.8mm Al
Power supply method	Battery	Battery
Remote exposure interface	RJ11	RJ11
image processing	none	Optional
Working environment temperature	-10°C~40°C	-10°C~40°C
size	221153.5261mm	225155263mm
weight	3.5kg	3.8kg